

Comprehensive Technical Analysis of Gambro AK 200 S and AK 200 ULTRA S Dialysis Systems

Published: January 23, 2026 | Index ID: #381

The evolution of renal replacement therapy (RRT) has been punctuated by the development of sophisticated delivery systems designed to mimic the intricate filtration processes of the human kidney. Among the most iconic and enduring platforms in this field are the **Gambro AK 200 S** and its advanced counterpart, the **AK 200 ULTRA S**. These machines represent a pinnacle of Swedish engineering, integrating complex fluid dynamics, electrochemical monitoring, and real-time software controls to provide life-sustaining treatment for patients with end-stage renal disease (ESRD).

As medical technology transitions toward more intensive therapies like Online Hemodiafiltration (HDF), understanding the underlying architecture of these devices becomes paramount for clinical engineers, nephrologists, and healthcare technicians. This guide provides an exhaustive technical breakdown of the AK 200 series, focusing on its operational mechanics, maintenance protocols, and the theoretical frameworks that govern its performance.

Theoretical Framework of Extracorporeal Therapy

To understand the AK 200 series, one must first grasp the physical principles it facilitates: **diffusion**, **convection**, and **ultrafiltration**. While the standard AK 200 S is primarily focused on high-efficiency hemodialysis (HD), the ULTRA S variant is specifically engineered to handle high-volume convection through Hemodiafiltration (HDF) and Hemofiltration (HF).

1. Diffusion (Solute Clearance)

Diffusion is the movement of solutes from an area of high concentration (blood) to low concentration (dialysate) across a semi-permeable membrane. The AK 200 S optimizes this by maintaining a counter-current flow, where blood and dialysate move in opposite directions to maximize the concentration gradient. The machine's ability to stabilize dialysate flow rates (Q_d) up to 1000 ml/min (though typically operated at 500 ml/min) is critical for effective urea and creatinine clearance.

2. Convection and Solvent Drag

In Hemodiafiltration, solutes are removed not just by concentration gradients but by 'dragging' them along with the fluid being removed. This is particularly effective for 'middle molecules' like β_2 -microglobulin. The **AK 200 ULTRA S** utilizes an integrated system to produce large volumes of sterile substitution fluid 'online' from the dialysate itself, avoiding the need for expensive, pre-packaged bags.

Technical Architecture: The Hydraulic and Electronic Core

The Gambro AK 200 S operates on a sophisticated hydraulic circuit that manages water purification, concentrate proportioning, and waste removal. The system is divided into several modules, each with dedicated sensors and actuators.

The Fluid Preparation System

The process begins with treated water entering the machine through an inlet valve. The water passes through a

Heat Exchanger, where it is pre-heated by the warm, used dialysate exiting the system, maximizing energy efficiency. It is then heated to the target temperature (typically 35–39°C) by a primary heater. A **Degassing Pump** creates a partial vacuum to remove dissolved air, preventing air bubbles from interfering with the sensors or entering the dialyzer.

Proportioning and Conductivity Control

The AK 200 S uses a volumetric proportioning system. It draws Acid and Bicarbonate concentrates using precise ceramic pumps. The resulting mixture's conductivity is measured by redundant sensors. If the conductivity deviates from the set window (usually 12.8 to 14.5 mS/cm), the **Bypass Valve** immediately diverts the fluid to the drain to prevent patient injury. This 'closed-loop' feedback system is the cornerstone of patient safety in the AK 200 series.

The Blood Circuit

The extracorporeal circuit consists of the arterial line, the blood pump, the dialyzer, and the venous line. The blood pump is a peristaltic roller pump designed to minimize hemolysis while maintaining precise flow rates (Q_b). Pressure transducers monitor **Arterial Pressure** (pre-pump), **Venous Pressure** (post-dialyzer), and **Transmembrane Pressure (TMP)**. TMP is a calculated value: $TMP = \frac{(P_{bi} + P_{bo})}{2} - \frac{(P_{di} + P_{do})}{2}$ where P represents pressure at various inlets (i) and outlets (o) for blood (b) and dialysate (d).

Comparative Analysis: AK 200 S vs. AK 200 ULTRA S

While sharing a common DNA, the ULTRA S model includes specialized components for advanced fluid management. The following table highlights the technical differences between the two models.

Feature	AK 200 S (Standard)	AK 200 ULTRA S
Primary Treatment Modes	HD, Isolated UF	HD, HDF (Online), HF (Online)
Substitution Fluid Source	External Bags (if HDF)	Online Production (U9000 Filters)
Fluid Balancing System	Volumetric (Balance Chambers)	Advanced Volumetric + Substitution Pump
Filtration System	Optional 1-Filter	Standard 3-Filter System (Ultra-pure)
Kt/V Monitoring	DIASCAN (Optional)	DIASCAN (Integrated)
Target Application	Routine Chronic Care	High-Efficiency/Critical Care

The 3-Filter Hemodialysis System

A defining characteristic of the **AK 200 ULTRA S** is its triple-filtration pathway. To ensure the safety of 'online' substitution fluid, the machine utilizes **Gambro U9000 filters**. These filters contain polyamide membranes that act as an absolute barrier to endotoxins and bacteria. The first filter purifies the incoming water, the second purifies the dialysate, and the third provides a final 'polishing' step for the substitution fluid before it enters the blood circuit. This ensures that the fluid meets the highest 'ultrapure' standards, significantly reducing systemic inflammation in chronic patients.

Maintenance and Service Protocols

According to the **Gambro AK 200 S Service Manual**, preventive maintenance is not merely a recommendation but a requirement for the device's certification. A qualified technician must perform a full calibration and safety check at least once per year or every 2,000 operating hours.

Step-by-Step Fluid Check (Fch) Procedure

The **Fch (Fluid Check)** is a diagnostic routine performed during the machine's power-on sequence. If the machine fails Fch, it will not enter the 'Ready' state.

1. Initial Power-On: Restart the machine and wait for the software (e.g., Program Version 8.xx) to initialize.
2. Fch Activation: The machine automatically cycles through various valve tests and pressure leak tests.
3. Manual Flow Set: As specified in the technical instructions, press the Fluid Flow button and set the main flow rate to 400 ml/min to verify the pump's calibration against the internal flow meter.
4. Calibration of Conductivity Cells: Technicians must use a reference meter to ensure the internal cells are within ± 0.1 mS/cm of the actual value.

Maintenance Checklist for Technicians

- Battery Replacement: The internal backup battery must be tested to ensure the machine can maintain the blood pump for at least 15-20 minutes during a power failure.
- O-Ring Inspections: Check all connector O-rings for signs of wear or calcification.
- Descaling and Disinfection: Regular chemical disinfection (using Peracetic Acid or Sodium Hypochlorite) and heat disinfection are vital to prevent biofilm formation in the hydraulic paths.
- Transducer Protector Check: Ensure that the air-fluid barriers in the pressure ports are dry and intact to prevent cross-contamination.

Troubleshooting and Operational Challenges

Operating a complex device like the AK 200 S involves managing various alarms. Understanding the root cause of these alarms is essential for reducing downtime.

Case Study: Low Conductivity Alarm

Observation: The machine enters 'Bypass' mode with a 'Cond Low' alarm.

Technical Analysis: This can be caused by an empty concentrate container, a kinked suction tube, or a malfunctioning concentrate pump. However, a common 'ghost' cause is a failure in the **Degassing System**. If air is present in the fluid path, the conductivity sensors will read incorrectly.

Solution: Check the vacuum level of the degassing pump. If the vacuum is too low, the air will not be removed, causing fluctuating conductivity readings. Adjust the degassing orifice or replace the pump motor if necessary.

Case Study: Blood Return Failure

Observation: During treatment, a power failure or technical error occurs, and the machine stops.

Technical Action: As per the 2015 safety advisories, users are instructed to **manually return blood** to the patient if the machine cannot be restarted immediately. This is done using the manual hand-crank on the blood pump. Technicians must ensure that the hand-crank is always accessible and that staff are trained in 'emergency manual return' to prevent blood loss in the circuit.

The Role of Software in Therapy Optimization

Modern versions of the AK 200 series (Version 8.xx and above) incorporate the **DIASCAN** system. This is a non-invasive tool that measures **Ionic Dialysance**. By slightly varying the conductivity of the dialysate and measuring the response, the machine can calculate the real-time clearance (Kt/V) and the total delivered dose ($\text{Kt/V}\cdot\text{h}$).

The Kt/V Calculation

The machine uses the following formula to track treatment adequacy:

By monitoring this continuously, the AK 200 S allows clinicians to adjust flow rates or treatment times mid-session to ensure the patient receives the prescribed dose, a process known as **Targeted Dialysis**.

Engineering for Safety and Longevity

The AK 200 series is renowned for its robust construction. The use of stainless steel in the heater vessels and ceramic in the proportioning pumps ensures that the machine can withstand thousands of hours of exposure to harsh chemicals and high temperatures. Furthermore, the **dual-processor architecture** provides a 'fail-safe' environment. One processor handles the control of the machine, while the second independently monitors all safety parameters. If the two processors disagree, the machine automatically enters a safe 'alarm' state.

Future Implications for Renal Technology

While the AK 200 series has been succeeded by newer platforms like the Artis or the AK 98, its legacy remains in the fundamental principles it perfected. The shift toward **High-Volume Online HDF**, which the ULTRA S pioneered, is now the standard of care in many parts of the world due to its superior cardiovascular outcomes. The rigorous maintenance standards and technical workflows established for the AK 200 continue to serve as the blueprint for training clinical engineers globally.

In conclusion, the Gambro AK 200 S and AK 200 ULTRA S are more than just dialysis machines; they are sophisticated medical instruments that require a deep understanding of both fluid mechanics and clinical physiology. By adhering to the strict maintenance protocols outlined in the service manuals and utilizing the advanced features like online filtration and DIASCAN monitoring, healthcare providers can ensure the highest level of safety and efficacy for their patients. The durability and precision of these systems continue to make them a benchmark in the field of extracorporeal therapy, proving that sound engineering and technical transparency are the keys to successful long-term patient outcomes.

Baca Juga (Artikel Terkait)

- **[Biomaterials Science: A Comprehensive Technical Analysis of Materials in Medicine](http://123caramba.com/biomaterials-science-technical-analysis-materials-medicine.pdf)**

URL: <http://123caramba.com/biomaterials-science-technical-analysis-materials-medicine.pdf>

- **[Comprehensive Guide to Biomedical Engineering: Essential Literatures, Core Technical Domains, and Practical Implementation](http://123caramba.com/comprehensive-guide-biomedical-engineering-literatures-technical-domains.pdf)**

URL: <http://123caramba.com/comprehensive-guide-biomedical-engineering-literatures-technical-domains.pdf>

- **[Biomedical Engineering: A Comprehensive Guide to Bridging Medicine and Technology](http://123caramba.com/biomedical-engineering-bridging-medicine-technology-guide.pdf)**

URL: <http://123caramba.com/biomedical-engineering-bridging-medicine-technology-guide.pdf>

- **[The Definitive Guide to Biomedical Engineering: Fundamentals, Signals, and System Integration](http://123caramba.com/biomedical-engineering-fundamentals-signals-system-integration-guide.pdf)**

URL: <http://123caramba.com/biomedical-engineering-fundamentals-signals-system-integration-guide.pdf>

Tags: **#Gambro AK 200 S #AK 200 ULTRA S #Dialysis Machine Maintenance #Hemodiafiltration #Clinical Engineering**

